

**GE Biomed**

Medical System

# BM 80

## Biological Microscope



## Product Specification:

Optical System:	Universal infinity anti-fungus optical system
Viewing Head:	30° inclined, 360° rotatable, Interpapillary distance 48 - 75 mm
Nosepiece:	Quadruple ball bearing nosepiece with positive click stops
Eyepiece:	Wide field EWF 10X /18 mm
Objectives:	BOECO Infinite corrected Plan Achromatic Objectives 4x/0.10, 10x/0.25, 40x/0.65 (S), 100x/1.25 (S, Oil)
Stage:	Large sized double layer mechanical specimen stage. 132 x 142 mm, Moving range 75 x 40 mm,
Focusing:	Coaxial coarse and calibrated fine focus control, fine division 0,001 mm, Coarse stroke 37,7mm/rotation fine stroke 0,2mm/rotation, mowing Range 20 mm safety autofocus stop unit.
Condenser:	Abbe condenser, N.A. 1.25 bright field
Illumination:	1W LED lamp, light intensity adjustable. Life time minimum 20000 hours

**Email:** [info@ge-biomed.co.kr](mailto:info@ge-biomed.co.kr)

**Phone:** +82(0)1067606469

**Address:** Incheon city, Namdong gu Gansok 3 dong,  
Gansok L.H.1 danji, 105 dong 1301, South Korea