

GE-600H COMPRESSOR NEBULIZER



- Light Indications
- Blue Light: During Operation
- Orange Light: Battery is Low Power
- Orange Blinking: Battery is Charging
- Blue Light Blinking: Battery has been fully Charged

Mesh Nebulizer Accessories



Mesh-Mouth Piece



Adult Mask



Child Mask



Usb Cable

Specification

Method Of Operation:	Active Vibrating Mesh Technology
Power Supply:	Lithium Titanate Battery (DC 4.8V Rechargeable)
Power Consumption:	Approx.2.0W
Auto Off:	10 Mins
Vibrating Frequency:	Approx.107 kHz+/-10%
Nebulization Rate:	0.25 MI/Min
Particle Size:Approx.	(88.08+1.91)% within 0.5 to 5um range
MMAD:	2.25±0.35um
Capacity Of Medication Cup :	Max 12ml
Noise Level:	≤ 50dBA(1 Meter Distance)
Battery Chargeing Voltage And Current:	5.0V 1A~2A
Battery Chargeing Time:	Approx.90 Mins
Battery Non-Stop Usage Time:	Approx.100 Mins
Lithium Titanate Battery Cycle Life:	Cycle Life:0.5C Approx.7000 Times, 1C Approx.4000 Times
Dimension:	Approx.L40 X W40 X H90 mm
Weight:Approx.110g	
Operating Temperature And Humidity:	10~40,30~85%RH,800~1060hPa
Storage And Delivery Condition:	— 20~60,20~75%RH,800~1060hPa
Safety Level:	Type B Class II

Incheon city, Namdong gu Gansok 3 dong, Gansok L.H.1 danji,
105 dong 1301, South Korea
Phone Number: +82(0)1067606469, Email: info@ge-
biomed.co.kr, www.ge-biomed.co.kr

Dealer

Global Medical Engineering (BD) Ltd.
17/2, Topkhana Road, Dhaka-1000
Phone: +88 02 9564225
E-mail. info@gmebd.com