

**GEBIOMED**

Medical System

# GE-600F COMPRESSOR NEBULIZER



Compressor Nebulizer

Compressor Nebulizer Helps to treat respiratory disease  
(as bronchial catarrh/pneumonia/asthmatic,etc)

- Long-life compressor
- High nebulization rate
- Quality aerosol every time
- Cost-effective and efficient
- Suitable to use in the occasion of home&hospital

Specification	
Power	AC 230V 50Hz,120V 60Hz,220V 50Hz/60Hz,240V 50Hz,220V 60Hz, etc
PowerConsumption	180VA or 1.5A
Nebulizer medication Capacity	6ml
Partical Size	275%within 0.5 to 5um range
MMAD	3.0um
Average Nebulization Rate	0.2ml/min
Noise Level	Below 65dBA(DIN EN 13544-1)
Operation Pressure Range	9.5~20 Psi/65.5Kpa~137.9Kpa
Compressor Pressure Range	30~50 Psi/206Kpa~345Kpa
Liter Flow Range	6~10 lpm

Nebulizer Accessory	
SH-B type Standard Accessories	Nebulizer bottle(6ml),Air tube(1.6m).Mouth piece,Adult mask,Child Mask,Air filters(5pcs)
Optional	Nose piece

Optional Nebulizer Accessory	
SH-A type Accessories:	Nebulizer bottle(8ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs)
Optional	Nose piece
SH-C type Accessories:	Nebulizer bottle(8ml),Air tube(1.6m).Mouth piece,Adult mask.Child Mask,Air filters(5pcs)
Optional	Nose piece
SH-D type Accessories:	Nebulizer bottle(8ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs)
Optional	Nose piece
SH-E type Accessories:	Nebulizer bottle(7ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs)
Optional	Nose piece

Incheon city, Namdong gu Gansok 3 dong, Gansok L.H.1 danji, 105 dong 1301, South Korea  
 Phone Number: +82(0)1067606469, Email: info@ge-biomed.co.kr,  
[www.ge-biomed.co.kr](http://www.ge-biomed.co.kr)

## Dealer

Global Medical Engineering (BD) Ltd.  
 17/2, Topkhana Road, Dhaka-1000  
 Phone: +88 02 9564225