

**GEBIOMED**

Medical System

# GE-600K COMPRESSOR NEBULIZER



Compressor Nebulizer

Compressor Nebulizer Helps to treat respiratory disease  
(as bronchial catarrh/pneumonia/asthmatic,etc)

- Long-life compressor
- High nebulization rate
- Quality aerosol every time
- Cost-effective and efficient
- Suitable to use in the occasion of home&hospital

| Specification                 |   |
|-------------------------------|---|
| Power                         | AC 230V 50Hz,120V 60Hz,220V 50Hz/60Hz,240V 50Hz,220V 60Hz,etc |
| PowerConsumption              | 180VA or 1.5A   |
| Nebulizer medication Capacity | 6ml   |
| Partical Size                 | 275%within 0.5 to 5um range                                   |
| MMAD                          | 3.0um   |
| Average Nebulization Rate     | 0.2ml/min   |
| Noise Level                   | Below 65dBA(DIN EN 13544-1)                                   |
| Operation Pressure Range      | 9.5~20 Psi/65.5Kpa~137.9Kpa                                   |
| Compressor Pressure Range     | 30~50 Psi/206Kpa~345Kpa                                       |
| Liter Flow Range              | 6~10 lpm  |

| Nebulizer Accessory            |  |
|--------------------------------|--|
| SH-B type Standard Accessories | Nebulizer bottle(6ml),Air tube(1.6m).Mouth piece,Adult mask,Child Mask,Air filters(5pcs) |
| Optional                       | Nose piece   |

| Optional Nebulizer Accessory |  |
|------------------------------|--|
| SH-A type Accessories:       | Nebulizer bottle(8ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs) |
| Optional                     | Nose piece   |
| SH-C type Accessories:       | Nebulizer bottle(8ml),Air tube(1.6m).Mouth piece,Adult mask.Child Mask,Air filters(5pcs) |
| Optional                     | Nose piece   |
| SH-D type Accessories:       | Nebulizer bottle(8ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs) |
| Optional                     | Nose piece   |
| SH-E type Accessories:       | Nebulizer bottle(7ml),Air tube(1.6m),Mouth piece,Adult mask,Child Mask,Air filters(5pcs) |
| Optional                     | Nose piece   |

Incheon city, Namdong guGansok 3 dong, Gansok L.H.1 danji,  
105 dong 1301, South Korea  
Phone Number: +82(0)1067606469, Email: info@ge-  
biomed.co.kr, www.ge-biomed.co.kr

## Dealer

Global Medical Engineering (BD) Ltd.  
17/2, Topkhana Road, Dhaka-1000  
Phone: +88 02 9564225  
E-mail. info@gmebd.com