

# Oxygen Concentrator

## GE W10L/B10L



Model Number: GE W10L/B10L

Operating Mode: button

Noise level 4dB

Power Consumption:400W

Oxygen flow: 0.5L/min~9.9L/min,Oxygen concentration 90%-96% (v/v)

Power supply: AC220V±10% 50Hz 60Hz(According to customer request)

Product size: 345(L)\*280(w)\*558(H)mm

Net weight 15.5kgs

Operating mode: button/remote control

LCD Display: present working time,total working time,presetting time from 1 min  
999mins.

Alarm: Cleaning filter alert.Built-in power failure alarm.compressor overheating  
stop alarm nebulizer(Atomization)



Nebulizer(Atomization )

Innovative Big

Flow Rate: 0.5-9.9 Lpm

Capacity(oxygen flow)	0.5L~9.9Lmin
Puity	93%±3%(5LPM)
Sound level	≤54db
Operation mode	Button
Nebulization function	
Net weight	15.5 kg
Size	345mmx280mmx558mm
Alarm	filter-clean reminder, built-in power failure alarm, Compressor overheating stop alarm
LED display	present working time,total working time.presetting time,filter cleaning remimder.fowrate data,purity data

Add:Incheon city, Namdong guGansok 3 dong, Gansok  
L.H.1 danji, 105 dong 1301, South Korea  
Phone: +82(0)1067606469  
info@ge-biomed.co.kr www.ge-biomed.co.kr

Dealer

Global Medical Engineering (BD) Ltd  
17/2,Topkhana Road, Dhaka-1000  
Bangladesh cell: +8801711020555  
info@gmebd.com www.gmebd.com